




**GLOBALWIN**

IMPLANT SYSTEM

# General Product Catalogue



**BIOSAFIN**

Partner of



Making Your Life Better.

# Table of Contents

## **GLOBALWIN Implant System**

Overview.....	p. 6- 7
GLOBALWIN Implants .....	p. 8-9
GLOBALWIN surgical kit .....	p. 10
Summary table for the use of GLOBALWIN burs.....	p. 11
Drills.....	p. 12-14
13mm Short Drill Stops.....	p. 15
Intra-operative Instruments and accessories .....	p. 15- 16
Surgical Kit and Prosthetics/Summary Table on Usage.....	p. 17 - 19
Prosthetic Accessories: Summary chart.....	p. 20

## **Prostheses with conventional procedure**

Implant closure: Implant cover screws .....	p. 21
Healing caps .....	p. 21

## **Impressions**

Repositioning Technique Impressions.....	p. 22
Pick-Up Technique Impressions .....	p. 23
Impression analogue .....	p. 24

## **Single cement/screw-retained prostheses**

Temporary abutments: temporary – rotating temporary .....	p. 24
Standard Abutments: Straight (0°) - Angled (15°) - Angled (25°) .....	p. 25
Aesthetic abutments: Straight (0°) - Angled (15°) - Angled (25°) .....	p. 26
Abutments for knife-edge technique .....	p. 27
Castable and Cast-on Abutments: .....	p. 28
– Castable in PMMA .....	p. 28
– Rotating castable in PMMA.....	p. 28
– for overcasting, in Platinor.....	p. 28
– rotating, for overcasting, in Platinor.....	p. 28
– for overcasting in CrCo .....	p. 28
– for overcasting in rotating CrCo.....	p. 28

## **Overdentures**

Ball attachment .....	p. 29
Analogue Ball attachment .....	p. 29
O-ring .....	p. 29
Box of O-Rings for ball attachment .....	p. 29
Analogue Ball attachment .....	p. 29

## **Multi-unit screw-retained prostheses**

Anti-rotational Trumpet Abutment with engagement and rotational Trumpet Abutment .....	p. 30
Performing Abutment: Straight (0°) - Angled (17°) - Angled (30°) - mounter .....	p. 31

## **Accessories for Trumpet and Performing Abutment**

Non-rotational Pick-Up Impression Coping for trumpet with engagement .....	p. 32
Performing Abutment rotational pick-up impression coping .....	p. 32
Healing cap for Performing Abutment .....	p. 32
Analogue for Performing Abutment and rotational and non-rotational trumpet.....	p. 32

## **Secondary Trumpet and Performing Abutments**

Anti-rotational Temporary Abutment for Trumpet.....	p. 33
Rotational Temporary Abutment for Performing Abutment and rotational Trumpet .....	p. 33
Anti-rotational Smooth Abutment for Performing Abutment.....	p. 33
Rotational Smooth Abutment for Performing Abutment.....	p. 33
Castable Abutment for Performing Abutment.....	p. 33

## **Prostheses with digital procedure**

Digital impression copings .....	p. 34
Digital impression copings for Performing Abutment .....	p. 34
Digital analogue .....	p. 34
Digital analogue for Performing Abutment .....	p. 34
Prosthetic Abutments: Base - Rotating - Preformed CAD CAM .....	p. 35

## **GLOBALWIN screws**

Screws and torque drivers .....	p. 37
---------------------------------	-------

## **GLOBALWIN SLIM Implant System**

Overview .....	p. 39- 40
Summary table for the usage of GLOBALWIN SLIM drills .....	p. 41
GLOBALWIN SLIM Implants.....	p. 42 - 43
SLIM surgical kit and prosthetics.....	p. 44

## **GLOBALWIN SLIM prostheses with conventional procedure**

Implant closure: Implant cover screws .....	p. 45
Healing caps .....	p. 45
Impressions .....	
Pick-Up Technique Impressions .....	p. 46
Impression analogue .....	p. 46
Single cement/screw-retained prostheses .....	
Temporary abutments: temporary – rotating temporary .....	p. 47
Standard Abutments: Straight (0°) - Angled (15°) - Angled (25°) .....	p. 48
Abutments for drilling.....	p. 49
Aesthetic abutments: Straight (0°) .....	p. 49
Aesthetic abutments: Angled (15°) - Angled (25°).....	p. 50
Abutments and cover screws for Flat Shift .....	p. 50
Castable rotating direct abutment.....	p. 51
Castable and Cast-on Abutments: .....	p. 51
– Rotating Flat Shift direct castable in PMMA .....	p. 51
– for overcasting in CrCo .....	p. 51

## **Overdentures**

Ball attachment .....	p. 52
Analogue Ball attachment .....	p. 52
O-ring.....	p. 52
Box of O-Rings for ball attachment .....	p. 52
Analogue Ball attachment .....	p. 52

## **Performing Abutments**

Globalwin SLIM straight (0°) and angled (17°) Performing Abutment .....	p. 53
Globalwin SLIM angled (30°) Performing Abutment / mounter .....	p. 54
Accessories for Performing Abutment and Trumpet .....	
Performing Abutment rotational pick-up impression coping .....	p. 55
Healing cap for Performing Abutment .....	p. 55
Analogue for Performing Abutment and rotational and non-rotational trumpet.....	p. 55

## **Secondary Performing and Trumpet Abutments**

Cement-retained rotational temporary abutment .....	
for Performing Abutment and rotational Trumpet.....	p. 56
Rotational Smooth Abutment for Performing Abutment.....	p. 56
Castable Abutment for Performing Abutment.....	p. 56

## **GLOBALWIN SLIM prostheses with digital procedure**

Digital impression copings for SLIM implants .....	p. 57
Digital impression copings for Performing Abutment .....	p. 57
Digital implant analogue .....	p. 57
Digital analogue for Performing Abutment .....	p. 57
Prosthetic Abutments: Anti-rotational base - Rotational - Preformed CAD CAM .....	p. 58

## **GLOBALWIN Guided surgery**

BIOSAFIN CREA Digital Centre .....	p. 59 - 61
Guided surgery kit instruments .....	p. 62
CBCT anatomical landmarks / Surgical Templates /	
Software and services / Accessories.....	p. 63

## **Restorative techniques for Multi-unit screw-retained prosthesis.....**

p. 65 - 67

## **Prosthesis Sterilisation Materials and Protocol .....**

p. 68

## **Patient Implant Passport .....**

p. 69

## **Select bibliography .....**

p. 70

# GLOBALWIN Implants

## The most versatile implants ever.



**GLOBALWIN** is the result of years of specialist experience from **BIOSAFIN**, an Italian company that in the past 15 years has taken the market of implant dentistry by storm, becoming a key player in the industry.

**GLOBALWIN** implants are the synthesis of clinical and laboratory expertise, concentrated into one single device that is ideal for the needs of today's dentist, combining the principles of ergonomics, surgical efficiency and prosthetic effectiveness into an intuitive and easy-to-use system.



### TRADEMARKS

- GLOBALWIN®
- MRS Micro Rough Surface®

### PATENTS REGISTERED

- Trumpet abutment: patent no. EP 3424460
- CAB – Clip Abutment Bar:
  - European Patent No 114 250327
  - International Patent PCT/EP2011/072448

**BIOSAFIN**

Partner of



Making Your Life Better.

# Global solution

*Micro and macro-morphology in perfect harmony*

*For all your prosthetics needs*



**Collar with micro-grooves**  
to lighten the stress at the crestal level.

**Gradual-depth spires**  
ensure maximum stability even with immediate-loading techniques.

**Truncated cone shape**  
with tapered apex, to facilitate the insertion of the implant and ensure maximum stability in different bone types


**MRS Micro Rough Surfaces**  
sandblasted and etched, the gold standard of osteoconductive surfaces. Clinical follow-up since 2006.

**Sharp apex**  
for facilitated insertion of the implant.

**Universal compatibility with internal hex**  
compatible with all major certified implant systems on the market. Self-centring for optimum fastening.

**Prosthetic platform**  
available in all major CAD CAM libraries.


**Innovative Prosthetic solutions**  
adapted for digital processes.



# GLOBALWIN Implants

## Universal LINE

Each implant is equipped with a  
**Cover screw**



**Implant diameter**



*Implant cover screw included*

*CSC40*

**Implant length**

**CODE**

<b>6 mm</b>	-
<b>8 mm</b>	GU375008MRS
<b>10 mm</b>	GU375010MRS
<b>11.5 mm</b>	GU375115MRS
<b>13 mm</b>	GU375013MRS
<b>15 mm</b>	GU375015MRS
<b>18 mm</b>	GU375018MRS







**One single CONNECTION**

Same shape and diameter

Differentiated platform

Implant



 <p><b>Ø 4.30</b></p> <p><i>CSC40</i></p> <p>CODE</p>	 <p><b>Ø 5.00</b></p> <p><i>CSC40</i></p> <p>CODE</p>	 <p><b>Ø 6.00</b></p> <p><i>CSC60</i></p> <p>CODE</p>
GU430006MRS	GU500006MRS	GU600006MRS
GU430008MRS	GU500008MRS	GU600008MRS
GU430010MRS	GU500010MRS	GU600010MRS
GU430115MRS	GU500115MRS	GU600115MRS
GU430013MRS	GU500013MRS	GU600013MRS
GU430015MRS	GU500015MRS	-
GU430018MRS	-	-

# GLOBALWIN Surgical Kit



**Torque wrench**  
DTW  
DTW75

**Non-torque wrench**  
MTW

**Manual adapter for non-torque wrench instruments**  
ADMTR

**Manual adapter for contra-angle screwdriver**  
HASD

**Short drill STOPS**

**Pilot Drills**  
STD1L  
STD1S

**Lance-shaped drill**  
PRD

**Contra-Angle**  
DEX

**Contra-angle screwdrivers**  
HSD20  
HSD25  
HSD30

**Manual screwdrivers**  
DSD12  
DSD16  
DSD26

**Mounter for Ball attachment**  
HMBLT

**Mounter for Performing Abutment**  
HMPRA

**Manual screwdrivers**  
MTRS  
MTRM  
MTRL

**Contra-angle implant drivers**  
HMTRS  
HMTRM  
HMTRL

**Dual-diameter drills**  
STD2L STD2S  
STD3L STD3S  
STD4L STD4S  
STD5L STD5S  
STD6L STD6S  
STD7L STD7S  
STD8L STD8S  
STD9  
STD10

**Direction indicator and depth gauge**  
DIRS  
DIRL

**Cortical drills**  
BP35  
BP37  
BP43  
BP50  
BP60

PRD	STD1L	STD2L	STD3L	STD4L	STD5L	STD6L	STD7L	STD8L	DIRS	DIRL
	2.00	2.20/2.60	2.60/3.00	3.00/3.40	3.40/3.80	3.80/4.20	4.20/4.60	4.60/5.00	5.00/5.40	5.40/5.80
DEX	STD1S	STD2S	STD3S	STD4S	STD5S	STD6S	STD7S	STD8S	STD9	STD10
	HSD20	HSD25	HSD30		HMTRS	HMTRM	HMTRL	BP37	BP43	
	DSD12	DSD16	DSD26		MTRS	MTRM	MTRL	BP50	BP60	
	HMBLT	HMPRA								

# Summary table for the use of GLOBALWIN drills

## SOFT BONE

Nominal dimensions	Ø MAX Endosseous	Ø Apical Spire	Drill diameters	Last Drill code	Drill for use in particularly corticalised bone
3.75	3.90	2.40	2.20/2.60	STD2	-
4.30	4.30	2.80	2.20/2.60	STD2	-
5.00	5.00	3.20	2.60/3.00	STD3	-
6.00	6.00	4.80	3.00/3.40	STD4	(BP60)

## MEDIUM BONE

Nominal dimensions	Ø MAX Endosseous	Ø Apical Spire	Drill diameters	Last Drill code	Drill for use in particularly corticalised bone
3.75	3.90	2.40	2.60/3.00	STD3	(BP37)
4.30	4.30	2.80	3.00/3.40	STD4	(BP43)
5.00	5.00	3.20	3.80/4.20	STD6	(BP50)
6.00	6.00	4.80	4.60/5.00	STD8	(BP60)

## HARD BONE

Nominal dimensions	Ø MAX Endosseous	Ø Apical Spire	Drill diameters	Last Drill code	Drill for use in particularly corticalised bone
3.75	3.90	2.40	3.00/3.40	STD4	BP37
4.30	4.30	2.80	3.40/3.80	STD5	BP43
5.00	5.00	3.20	4.20/4.60	STD7	BP50
6.00	6.00	4.80	5.00/5.40	STD9	BP60


# Drills

Made of:


*in surgical steel, sterilisable*

## Drill features:

- dual diameter
- variable length



## Ratio between: Drill - Implant - Direction indicator and depth gauge



The drill lengths are overestimated by 0.5 mm compared to the endosseous length of the implants.

## Recommendations:

- For reasons of safety, it is advisable to replace the drills after a maximum of 50 cuts.
- In any case, check the cutting efficiency before each use, taking into account that in compact bone there is greater wear.


# Drills

Made of:

*in surgical steel, sterilisable*

## Lance-shaped drill

Item code
PRD




## Pilot Drill

Diameter	Length	Item code
Ø 2.00 mm	6-13 mm	STD1S
Ø 2.00 mm	5-18 mm	STD1L



## Dual-diameter long drill

Diameter	Length	Item code
Ø 2.20/2.60 mm	5-18 mm	STD2L
Ø 2.60/3.00 mm	5-18 mm	STD3L
Ø 3.00/3.40 mm	5-18 mm	STD4L
Ø 3.40/3.80 mm	5-18 mm	STD5L
Ø 3.80/4.20 mm	5-18 mm	STD6L
Ø 4.20/4.60 mm	5-18 mm	STD7L
Ø 4.60/5.00 mm	5-18 mm	STD8L




# Drills

Made of:

***in surgical steel, sterilisable***


## Dual-diameter short drill

Diameter	Length	Item code
Ø 2.20/2.60 mm	6–13 mm	STD2S
Ø 2.60/3.00 mm	6–13 mm	STD3S
Ø 3.00/3.40 mm	6–13 mm	STD4S
Ø 3.40/3.80 mm	6–13 mm	STD5S
Ø 3.80/4.20 mm	6–13 mm	STD6S
Ø 4.20/4.60 mm	6–13 mm	STD7S
Ø 4.60/5.00 mm	6–13 mm	STD8S
Ø 5.00/5.40 mm	6–13 mm	STD9
Ø 5.40/5.80 mm	6–13 mm	STD10



## Cortical Drill

Diameter	Item code
Ø 3.40 mm	BP35
Ø 3.65 mm	BP37
Ø 4.20 mm	BP43
Ø 4.90 mm	BP50
Ø 5.90 mm	BP60



*IT is necessary to use the bone-finishing instrument with the corresponding 5 mm-high healing cap.*

# 13mm short drill stop

Made of:

***in surgical steel, sterilisable***

## Short DRILL STOPS

STD1S	STG6N
STD2S	STG8N
STD3S	STG10N
	STG115N
STD4S	STG6R
STD5S	STG8R
STD6S	STG10R
	STG115R
STD7S	STG6W
STD8S	STG8W
STD9	STG10W
STD10	STG115W



# Intra-operative Instruments and Accessories

Made of:


***Gr5 titanium, sterilisable***

## Direction indicator and depth gauge

	Item code
Short	DIRS
Long	DIRL

**USE:** after Ø 2.0 mm drill


**TO:** check the depth of the surgical site



## Surgical guide

*to control the inclination of the drill from 0° to 45°.  
Can be completely dismantled for washing and sterilisation.*

	Item code
	SG46



# Intra-operative Instruments and Accessories

Made of:

***in surgical steel, sterilisable***

## Visual Card

Item code
VCGU

**USE:** during treatment design, to be superimposed onto patient x-rays

**TO:** identify the size of the implant to use

**DIMENSIONS TO SCALE:** 1:1 - 1.1:1 - 1.2:2

## GLOBALWIN autoclavable TRAY

	Item code
Standard Plus	GWST
Standard	GWSTS




# Surgical and prosthetic instruments

Made of:

***titanium, sterilise before use***

## Non-torque surgical wrench

Item code
MTW






# Surgical and prosthetic instruments

Made of:

***in surgical steel, sterilisable***

## Manual implant driver

	Item code
Short	MTRS
Medium	MTRM
Long	MTRL



## Contra-angle implant driver

	Item code
Short	HMTRS
Medium	HMTRM
Long	HMTRL



## Manual adapter for non-torque wrench instruments

(codes MTRS and MTRL)

	Item code
	ADMTR




# Surgical and prosthetic instruments

Made of:

***in surgical steel, sterilisable***

## Torque wrench

	Item code
Torque 10-35 Ncm	DTW
Torque 10-75 Ncm	DTW75



## DTW75 torque wrench adaptor for MTR series manual instruments

	Item code
	DTWM



## Contra-angle screwdriver

*Also suitable for DTW and HASD*

Length	Item code
20 mm	HSD20
25 mm	HSD25
30 mm	HSD30



## Manual adapter for contra-angle instruments

	Item code
	HASD



# Surgical and prosthetic instruments

Made of:

***in surgical steel, sterilisable***

## Contra-Angle


	Item code
	DEX



## Manual screwdriver

*Also suitable for use with torque wrench*

Length	Item code
12 mm	DSD12
16 mm	DSD16
26 mm	DSD26



# Summary table of usage

<b>DTW</b>	+	→ HASD	+	HSD20	* If necessary, use with the DEX extension
		→ DSD12 → DSD16 → DSD26		HSD25 HSD30	
<b>DTW 75</b>	→	DTWM	+	99SC	
		HASD		99MC	
		DSD12 DSD16 DSD26		99LC	
<b>MTW</b>	→	99SC		HSD20	
		99MC		HSD25	
		99LC		HSD30	

# Prosthetics Accessories - Summary chart

Temporary Prostheses		TTA40	TTA40R													
Cement-retained prostheses		SAF	ASA	AAA 15°	AAA 25°	SSA	SAA 15°	SAA 25°	OCA	CSA	CrA	BSA	Custom	PMA		
Screw-retained prostheses Single		SAF	ASA	SSA			OCA			CSA	CrA	BSA	Custom	PMA		
Screw-retained prostheses Multiple	onto PRA abutments	PRA 0°	PRA 17°	PRA 30°	directly onto Implants										Welding	directly onto Implants
CAB	onto PRA abutments	PRA 0°	PRA 17°	PRA 30°											TRU03	TRU03R
Welding accessories	onto PRA abutments	PRA 0°	PRA 17°	PRA 30°	PRAPA										Overdentures	BLT
		PRA 0°	PRA 17°	PRA 30°											TRU03	TRU03R

# Implant sealing:

Screws made from:



**titanium, sterilise before use**



## Cover screw



	Item code	Item code
	CSC40	CSC60



## Healing cap



Gingival height	Item code	Item code	Item code	Item code
GH 3.0 mm	HSC4003	HSC5003	HSC6003	HSC7503
GH 5.0 mm	HSC4005	HSC5005	HSC6005	HSC7505
GH 7.0 mm	HSC4007	HSC5007	HSC6007	HSC7507


# Repositioning Technique Impressions

ICSSC screw for  
repositioning technique  
impression coping



Abutments made from:

**titanium, sterilise before use**



## Impression coping

ICSSC prosthetic screw included




	Item code	Item code	Item code	Item code
	ICS40	ICS50	ICS60	ICS75

## Rotational impression coping

ICSSC prosthetic screw included

Item code
ICSR





# Pick-Up Technique Impressions

ICPSC screw for pick-up impression coping



Abutments made from:

**titanium, sterilise before use**




## Impression coping

ICPSC prosthetic screw included

Item code	Item code	Item code
ICP40	ICP50	ICP60

## Impression coping for pressure technique

Item code
ICPF



# Impression analogue

Made of:  
***in titanium***

## Analogue

	<b>Item code</b>
	ANG



# Temporary Abutments

Abutments made from:  
***titanium, sterilise before use***


*PSC screw for  
prosthetic  
abutment*



## Temporary Abutment

*PSC prosthetic screw included*


	<b>Item code</b>
	TTA40



## Rotational Temporary Abutment

*PSC prosthetic screw included*

	<b>Item code</b>
	TTA40R





# Standard Abutments

Abutments made from:

**titanium, sterilise before use**


PSC screw for  
prosthetic  
abutment



## 0° Straight Standard Abutment

*PSC prosthetic screw included*


Diameter	Gingival height	Item code
Ø 3.75 mm	GH 2 mm	SSA3750020
Ø 3.75 mm	GH 4 mm	SSA3750040
Ø 4.50 mm	GH 2 mm	SSA450020
Ø 4.50 mm	GH 4 mm	SSA450040
Ø 5.00 mm	GH 2 mm	SSA500020
Ø 5.00 mm	GH 4 mm	SSA500040
Ø 6.00 mm	GH 2 mm	SSA600020
Ø 6.00 mm	GH 4 mm	SSA600040



## 15° Angled Standard abutment

*PSC prosthetic screw included*


Diameter	Gingival height	Item code
Ø 5.00 mm	GH 2.00 mm	SAA501520
Ø 5.00 mm	GH 4.00 mm	SAA501540



## 25° Angled Standard abutment

*PSC prosthetic screw included*

Diameter	Gingival height	Item code
Ø 5.00 mm	GH 2.00 mm	SAA502520
Ø 5.00 mm	GH 4.00 mm	SAA502540



# Aesthetic abutments:

Abutments made from:

**titanium, sterilise before use**


PRAPSC screw  
for prosthetic  
abutment



## 0° Straight Aesthetic Abutment

PRAPSC prosthetic screw included


Diameter	Gingival height	Item code
Ø 3.75 mm	GH 2 mm	ASA3750020
Ø 3.75 mm	GH 4 mm	ASA3750040
Ø 4.50 mm	GH 2 mm	ASA450020
Ø 4.50 mm	GH 4 mm	ASA450040
Ø 5.00 mm	GH 2 mm	ASA500020
Ø 5.00 mm	GH 4 mm	ASA500040
Ø 6.00 mm	GH 2 mm	ASA600020
Ø 6.00 mm	GH 4 mm	ASA600040



## 15° Angled Aesthetic Abutments

PRAPSC prosthetic screw included


Diameter	Gingival height	Item code
Ø 5.00 mm	GH 1.50 mm	AAA501515
Ø 5.00 mm	GH 3.00 mm	AAA501530



## 25° Angled Aesthetic Abutments

PRAPSC prosthetic screw included

Diameter	Gingival height	Item code
Ø 5.00 mm	GH 1.50 mm	AAA502515
Ø 5.00 mm	GH 3.00 mm	AAA502530



# Abutments for knife-edge technique

PSC screw for  
prosthetic  
abutment




Abutments made from:

**titanium, sterilise before use**

## Straight abutment for knife-edge technique

*PSC prosthetic screw included*


Diameter	Item code
Ø 4.00 mm	SAF40
Ø 5.00 mm	SAF50
Ø 6.00 mm	SAF60



## Straight abutment for rotating knife-edge technique

*PSC prosthetic screw included*

	Item code
Ø 4.00 mm	SAF40R
Ø 5.00 mm	SAF50R



# Castable and Cast-On Abutments

PSC screw for  
prosthetic  
abutment




**sterilise before use**

## Castable PMMA Abutment

*PSC prosthetic screw included*


Item code
CSA40



## Castable Rotational PMMA Abutment

*PSC prosthetic screw included*


Item code
CSA40R



## Platinor Overcasting Abutment

*PSC prosthetic screw included*


Item code
OCA40



## Platinor Rotational Overcasting Abutment

*PSC prosthetic screw included*


Item code
OCA40R



## CrCo Overcasting Abutment

*PSC prosthetic screw included*


Item code
CrA40



## CrCo Rotational Overcasting Abutment

*PSC prosthetic screw included*

Item code
CrA40R




# Overdentures

Attachments, analogues, o-rings and boxes made from:  
**titanium, sterilise before use**

## Ball attachment

Box of O-Rings for ball attachment BLTBOX included

Gingival height	Item code
GH 1.5 mm	BLT4015
GH 3.0 mm	BLT4030
GH 5.0 mm	BLT4050



## Analogue Ball attachment

Item code
BLTA



## O-ring

Item code
BLTOR



## Box of O-Rings for ball attachment

Item code
BLTBOX




## Mounter for Ball attachment

Item code
HMBLT



# Trumpet Abutment multi-unit screw-retained prosthesis



Abutments made from:

**titanium, sterilise before use**


Trumpet patented abutment:  
patent no. EP 3424460

**THE PATENTED TRUMPET ABUTMENT** expands the available surface of the soft tissues, promoting the creation of an ideal bacterial seal. It's specifically designed to ensure a perfect seal, thanks to the articulated attachment of soft tissue to the neck of the abutment.

## Anti-rotational Trumpet abutment with engagement


*TRUPSC prosthetic screw included*

Item code
TRU03



## Rotational Trumpet Abutment

Item code
TRU03R



## Prosthetic screw for Trumpet Abutment

Item code
TRUPSC



# Performing Abutment multi-unit screw-retained prosthesis

PRAPSC screw for secondary abutment



Abutments made from:

**titanium, sterilise before use**




Pre-assembled Mounter included in the pack of straight Performing Abutments for correct placement and initial tightening turns.



Tightening of the Abutment must then be carried out at  $\leq 30$  N/cm with the Mounter code HMPRA

## Performing Abutment:

Angulation	Gingival height	Item code
00°	GH 1.5 mm	PRA0015
00°	GH 2.5 mm	PRA0025
00°	GH 3.0 mm	PRA0030




PRAPSC prosthetic screw included

17°	GH 2.5 mm	PRA1725
17°	GH 3.5 mm	PRA1735



PRAPSC prosthetic screw included

30°	GH 2.5 mm	PRA3025
30°	GH 3.5 mm	PRA3035



## Mounter for Performing Abutment

Item code
HMPRA



Pre-assembled mounter included in the pack of angled abutments for placement and verification of the angle.



# Accessories for Trumpet and Performing Abutment


Abutments made from:

**titanium, sterilise before use**

## Non-rotational Pick-Up Impression Coping for Trumpet with engagement

*Screw for PRAICPSC secondary abutment included*


Item code
PRAICP



## Performing Abutment rotational Pick-up impression coping

*Screw for PRAICPSC secondary abutment included*

Item code
PRAICPR



## Healing cap for Performing Abutment

Gingival height	Item code
GH 3.50 mm	PRAHSC35
GH 6.00 mm	PRAHSC60



## Analogue for Performing Abutment and rotational and non-rotational Trumpet

Item code
PRAA





# Secondary Trumpet and Performing Abutments

PRAPSCS Screw




Abutments made from:

**titanium, sterilise before use**

## Anti-rotational Temporary Abutment for Trumpet

*Screw for PRAPSCS secondary abutment included*


Item code
PRATA



## Rotational Temporary Abutment and/or cementing link for Performing Abutment and rotational Trumpet

*Screw for PRAPSCS secondary abutment included*

Item code
PRATAR*



## Performing Abutment Anti-rotational Smooth Abutment

*Screw for PRAPSCS secondary abutment included*

Item code
PRAPA



## Performing Abutment Rotational Smooth Abutment

*Screw for PRAPSCS secondary abutment included*


Item code
PRAPAR



## Performing Abutment Castable Abutment

*Screw for PRAPSCS secondary abutment included*

Item code
PRACS



\*can also be used for CAD-CAM modelling by removing the retentive part

# Prostheses with digital procedure

Abutments made from:


**titanium, sterilise before use**

## Digital impression coping \*\*

*PSC prosthetic screw included*

*It is advisable to tighten the supplied screw to 10 Ncm.*

Item code
ICD*




## Digital impression coping \*\* for Performing Abutment

*Screw for PRAPSCS secondary abutment included*

*It is advisable to tighten the supplied screw to 10 Ncm.*


Item code
ICDPRA*



\*autoclavable


## Digital implant analogue \*\*

Item code
DAI



## Digital analogue for Performing Abutment

Item code
DAPRA



*\*can also be used for CAD-CAM modelling by removing the retentive part*

*\*\*Libraries can be downloaded from [www.globalwin.eu](http://www.globalwin.eu)*

# Prosthetic abutments

Abutments made from:

**titanium, sterilise before use**

PSC screw for  
prosthetic  
abutment



## Anti-rotational base abutment \*\*

*PSC prosthetic screw included*

Gingival height	Item code
GH 1.5 mm	BSA4015
GH 3.0 mm	BSA4030



## Rotational Base Abutment \*\*

*PSC prosthetic screw included*


Gingival height	Item code
GH 1.5 mm	BSA4015R
GH 3.0 mm	BSA4030R



## Premilled Abutments\*\*

*PSC prosthetic screw included*

Diameter	Item code
Ø 10.00 mm	PMA10
Ø 14.00 mm	PMA14
In cobalt-chrome	
Ø 10.00 mm	PMA10CR
Ø 14.00 mm	PMA14CR


















\*\*Libraries can be downloaded from [www.globalwin.eu](http://www.globalwin.eu)



**GLOBALWIN**  
IMPLANT SYSTEM


# GLOBALWIN screws

Description	Item code	Tightening torque	
Cover Screw	<b>CSC</b>	10 Ncm	
Healing cap	<b>HSC</b>	10 Ncm	
Screw for Standard impression coping	<b>ICSSC</b>	10 Ncm	
Screw for pick-up impression coping	<b>ICPSC</b>	10 Ncm	
Prosthetic Screw	<b>PSC</b>	30 Ncm	
Performing Abutment angled Prosthetic Screw	<b>PRAPSC</b>	30 Ncm	
Secondary prosthetic screw	<b>PRAPSCS</b>	20 Ncm	
Healing cap for Performing Abutment	<b>PRAHSC</b>	10 Ncm	
Prosthetic Screw for Trumpet Abutment	<b>TRUPSC</b>	30 Ncm	
Screw for impression coping	<b>PRAICPSC</b>	10 Ncm	
Digital analogue screws	<b>WDAV</b>	20 Ncm	
Slim prosthetic screw	<b>VP2SI</b>	30 Ncm	
Screw for Slim impression coping	<b>VP3SI</b>	30 Ncm	
Prosthetic screw for Flat Shift	<b>VFS</b>	25 Ncm	
Screws for transfer	<b>VMILSI</b>	10 Ncm	




**GLOBALWIN**  
IMPLANT SYSTEM

# SLIM IMPLANTS 2.9 AND 3.3





# GLOBALWIN SLIM INDICATIONS FOR INDICATIONS FOR IMPLANTS Ø 2.9 AND Ø 3.3




- for central and lateral mandibular incisors and maxillary lateral incisors
- for very limited interproximal spaces
- to avoid bone regeneration techniques
- in case of orthodontic therapies
- in case of thin crests

**A: Static Structural**  
Equivalent Stress  
Type: Equivalent (von-Mises) Stress  
Unit: MPa  
Time: 1



Worst Case test carried out on a SLIM BIOSAFIN implant.



The test, carried out by simulating the WORST CASE of a BIOSAFIN SLIM 2.9 implant with the GH3 straight abutment (a condition that does not occur in clinical practice, thanks to the use of angled abutments), highlights the Von Mises stress performance, with a load of 500N, obtained with the Finite Element method.

The system amply meets the chewing load levels for which it was designed.



## SLIM IMPLANT Ø 2.9 AND Ø 3.3

GLOBALWIN SLIM 2.9 and 3.3 implants with SLIM internal hex connection extends and completes the GLOBALWIN Line of Implants, to meet the needs of Dentists for implant placement in very narrow anatomical spaces. These implants are only indicated for the mandibular central and lateral incisors and maxillary lateral incisors.

### Micro- and macro-morphology in perfect harmony

**Neck with micro-grooves** to lighten the stress at the crestal level.

**Gradual-depth threads** ensure maximum stability even with immediate-loading techniques.

**Truncated cone shape** with tapered apex, to facilitate insertion of the implant and ensure maximum stability in different bone types.

**MRS Micro Rough Surface** sandblasted and etched, the gold standard of osteoconductive surfaces. Clinical follow-up since 2005.


**Sharp apex** for facilitated insertion of the implant.

### SLIM internal hex connection with Free Lock technology

**Innovative Prosthetic solutions** adapted for digital processes.



- Used since 1995.
- Helps achieve great connection stability



GLOBALWIN SLIM implants are one of the smallest on the market and thanks to the internal connection with Free Lock technology, used since 1995, they provide excellent stability and load resistance.

Its cylindrical shape, with a self-tapping conical apex, MRS surface and microgrooved neck make it a high-performance implant.

- Ø 1.40 mm in implant 2.9
- Ø 1.50 mm in implant 3.3

## FULL RANGE OF

### SLIM implants

**SLIM CONNECTION**



- 29010GSI/MRS
- 29115GSI/MRS
- 29013GSI/MRS
- 29015GSI/MRS



- 33010GSI/MRS
- 33115GSI/MRS
- 33013GSI/MRS
- 33015GSI/MRS

### DIAMETERS AND HEIGHTS

Ø 2.9	Ø 3.3
h 10	h 10
h 11.5	h 11.5
h 13	h 13
h 15	h 15



# Summary table for the usage of GLOBALWIN SLIM drills

## SOFT BONE

Nominal dimensions	Ø MAX Endosseous	Ø Apical Spire	Drill diameters	Last Drill code	Drill for use in particularly corticalised bone
2.90	3.10	1.80	2.20/2.60	STD2	-
3.30	3.50	1.80	2.20/2.60	STD2	-

## MEDIUM BONE

Nominal dimensions	Ø MAX Endosseous	Ø Apical Spire	Drill diameters	Last Drill code	Drill for use in particularly corticalised bone
2.90	3.10	1.80	2.20/2.60	STD2	-
3.30	3.50	1.80	2.20/2.60	STD2	-


## HARD BONE

Nominal dimensions	Ø MAX Endosseous	Ø Apical Spire	Drill diameters	Last Drill code	Drill for use in particularly corticalised bone
2.90	3.10	1.80	2.20/2.60	STD2	-
3.30	3.50	1.80	2.20/2.60	STD2	BP35

Use instrument 62LSI to remove the SLIM implant from the ampoule.  
 After placement in the prepared site, proceed with screw-tightening with the dedicated moulder from the following instruments:  
 62SC; 62MC; 62LC; 99SC; 99MC; 99LC.

# GLOBALWIN SLIM IMPLANTS

Each implant is equipped with a  
**Cover screw**



## One single CONNECTION

Same shape and diameter

Differentiated platform

Implant

Ø 2.9 mm

Ø 3.3 mm



## Implant diameter



*Implant cover screw included*

VT29SI

VT33SI

Implant length

CODE

CODE

10 mm

29010GSI

33010GSI

11.5 mm

29115GSI

33115GSI

13 mm

29013GSI

33013GSI

15 mm

29015GSI

33015GSI

# GLOBALWIN SLIM Surgical kit and prosthetic instruments

Made of:

*in surgical steel, sterilisable*

## instrument for extraction of the implant from the ampoule

Item code
62LSI



## Manual implant driver

	Item code
Short	62SC
Medium	62MC
Long	62LC



## Contra-angle implant driver

	Item code
Short	99SC
Medium	99MC
Long	99LC



# Implant sealing:

Screws made from:


**titanium, sterilise before use**



## Cover screw



	Item code	Item code
	VT29SI	VT33SI



## Flared healing caps



Gingival height	Item code	Item code
GH 2.0 mm	VG292SI	VG332SI
GH 3.0 mm	VG293SI	VG333SI
GH 5.0 mm	VG295SI	VG335SI
GH 7.0 mm	VG297SI	VG337SI

# Pick-Up Technique Impressions

VMILSI screw for pick-up impression coping



Abutments made from:

**titanium, sterilise before use**



## Impression coping

VMILSI prosthetic screw included



	Item code	Item code
	IMG29LSI	IMG33LSI

# Impression analogue

Made of:

**in titanium**

## Analogue



	Item code	Item code
	A29SI	A33SI



# Temporary Abutments

Abutments made from:

**titanium, sterilise before use**


VP2SI screw for  
prosthetic  
abutment



## Temporary Abutment

VP2SI screw included


Gingival height	Item code	Item code
GH 1.0 mm	● MP2901SI	● MP3301SI
GH 2.0 mm	● MP2902SI	● MP3302SI
GH 3.0 mm	● MP2903SI	● MP3303SI



## Rotational Temporary Abutment (also suitable for welding)

VP2SI prosthetic screw included

Gingival height	Item code
GH 1.0 mm	● MDS331WRSI
GH 2.0 mm	● MDS332WRSI
GH 3.0 mm	● MDS333WRSI



# Standard Abutments

Abutments made from:

**titanium, sterilise before use**


VP2SI screw for  
prosthetic  
abutment



## Straight Standard abutment

VP2SI prosthetic screw included


Gingival height	Item code
GH 1.0 mm	● MD2901SI
GH 2.0 mm	● MD2902SI
GH 3.0 mm	● MD2903SI
GH 1.0 mm	● MD3301SI
GH 2.0 mm	● MD3302SI
GH 3.0 mm	● MD3303SI



## 15° Angled Standard abutment

VP2SI prosthetic screw included


Gingival height	Item code
GH 1.0 mm	● MIN29151SI
GH 2.0 mm	● MIN29152SI
GH 3.0 mm	● MIN29153SI
GH 1.0 mm	● MIN33151SI
GH 2.0 mm	● MIN33152SI
GH 3.0 mm	● MIN33153SI



## 25° Angled Standard abutment

VP2SI prosthetic screw included

Gingival height	Item code
GH 1.0 mm	● MIN33251SI
GH 2.0 mm	● MIN33252SI
GH 3.0 mm	● MIN33253SI





# Abutments for drilling

Abutments made from:

**titanium, sterilise before use**


## Straight Standard abutment

*VP2SI prosthetic screw included*

### Item code

● MF2900SI

● MF3300SI



# Aesthetic abutments:

Abutments made from:


**titanium, sterilise before use**

## Straight Aesthetic Abutment

*VP2SI prosthetic screw included*

Gingival height	Item code
GH 1.0 mm	● WE29001SI
GH 2.0 mm	● WE29002SI
GH 3.0 mm	● WE29003SI
GH 1.0 mm	● WE33001SI
GH 2.0 mm	● WE33002SI
GH 3.0 mm	● WE33003SI


*VP2SI screw for prosthetic abutment*



## 15° Angled Aesthetic Abutments

VP2SI prosthetic screw included


Gingival height	Item code
GH 1.0 mm	● WE29151SI
GH 2.0 mm	● WE29152SI
GH 3.0 mm	● WE29153SI
GH 1.0 mm	● WE33151SI
GH 2.0 mm	● WE33152SI
GH 3.0 mm	● WE33153SI



## 25° Angled Aesthetic Abutments

VP2SI prosthetic screw included

Gingival height	Item code
GH 1.0 mm	● WE33251SI
GH 2.0 mm	● WE33252SI
GH 3.0 mm	● WE33253SI



# Flat Shift abutments

Abutments made from:


**titanium, sterilise before use**

VFS screw for  
Flat Shift

## Flat Shift abutments

VFS prosthetic screw included

Gingival height	Item code
GH 3.0 mm	● MFS333SI
GH 4.0 mm	● MFS334SI
GH 5.0 mm	● MFS335SI



## Flat Shift cover screws

Item code

● VTFS33



# Cylindrical direct rotational abutments

Abutments made from:

**titanium, sterilise before use**


VP2SI screw for  
prosthetic  
abutment



## Flat Shift abutments

VP2SI prosthetic screw included

Gingival height	Item code
GH 1.0 mm	● MDC331SI
GH 2.0 mm	● MDC332SI
GH 3.0 mm	● MDC333SI



# Castable and Cast-On Abutments

**sterilise before use**


VP2SI screw for  
prosthetic  
abutment



## Rotating Flat Shift direct castable abutment in PMMA

VP2SI prosthetic screw included


Gingival height	Item code
GH 1.0 mm	MDF331SI



## CrCo Overcasting Abutment

VP2SI prosthetic screw included

Item code
● CR29SI
● CR33SI




# Overdentures

Attachments, analogues, o-rings and boxes made from:  
**titanium, sterilise before use**

## Ball attachment

Retentive O-ring and Box in titanium, Ø 2.1mm to be tightened with HMBLT instruments

Gingival height	Item code
GH 1.0 mm	MOR331SI
GH 2.0 mm	MOR332SI
GH 3.0 mm	MOR333SI
GH 4.0 mm	MOR334SI
GH 5.0 mm	MOR335SI



## Analogue Ball attachment

Item code
BLTA




## O-ring

Item code
BLTOR



## Box of O-Rings for ball attachment

Item code
BLTBOX



## Mounter for Ball attachment

Item code
HMBLT



# Performing Abutment:

Abutments made from:

**titanium, sterilise before use**



Pre-assembled Moulder included in the pack of straight Performing Abutments for correct placement and initial tightening turns.




Tightening of the Abutment must then be carried out at  $\leq 30$  N/cm with the Moulder code HMPRA

## Globalwin SLIM straight Performing Abutment

Pre-assembled moulder included, tighten to 30 N/cm

Gingival height	Item code
GH 1.0 mm	● EA33001SI
GH 2.0 mm	● EA33002SI
GH 3.0 mm	● EA33003SI
GH 4.0 mm	● EA33004SI
GH 5.0 mm	● EA33005SI



## Globalwin SLIM 17° angled Performing Abutment

Pre-assembled moulder included, tighten to 30 N/cm

Gingival height	Item code
GH 2.0 mm <i>VP3SI prosthetic screw included</i>	● EA33172SI
GH 3.0 mm <i>VP2SI prosthetic screw included</i>	● EA33173SI
GH 4.0 mm <i>VP2SI prosthetic screw included</i>	● EA33174SI



# Performing Abutment:


Abutments made from:

**titanium, sterilise before use**

## Globalwin SLIM 30° angled Performing Abutment

*Pre-assembled mounter included, 30 Ncm torque,  
VP2SI prosthetic screw included*

Gingival height	Item code
GH 4.0 mm	● EA33304SI
GH 5.0 mm	● EA33305SI



## Mounter for Performing Abutment

Item code
HMPRA



*Pre-assembled mounter included in the pack of angled abutments for  
placement and verification of the angle.*



# Performing Abutment Accessories


Abutments made from:

**titanium, sterilise before use**

## Performing Abutment rotational Pick-up impression coping

*Screw for PRAICPSC secondary abutment included*

Item code
PRAICPR



## Healing cap for Performing Abutment

Gingival height	Item code
GH 3.5 mm	PRAHSC35
GH 6.0 mm	PRAHSC60



## Analogue for Performing Abutment and rotational and non-rotational Trumpet

Item code
PRAA



# Performing Abutment and Trumpet secondary abutments

PRAPSCS Screw




Abutments made from:

**titanium, sterilise before use**

## Rotational Temporary Abutment and/or cementing link for Performing Abutment and rotational Trumpet

*Screw for PRAPSCS secondary abutment included*

Item code
PRATAR*



## Performing Abutment Rotational Smooth Abutment

*Screw for PRAPSCS secondary abutment included*


Item code
PRAPAR



## Performing Abutment Castable Abutment

*Screw for PRAPSCS secondary abutment included*

Item code
PRACS



*\*can also be used for CAD-CAM modelling by removing the retentive part*



# Prostheses with digital procedure

Abutments made from:

**titanium, sterilise before use**

## Digital impression copings for SLIM implants Ø2.9 and Ø3.3 \*\*

*PSC prosthetic screw included*

*It is advisable to tighten the supplied screw to 10 Ncm.*

Item code
● WS29HSI




## Digital impression coping for Performing Abutment

*Screw for PRAPSCS secondary abutment included*

*It is advisable to tighten the supplied screw to 10 Ncm.*


Item code
ICDPRA*



\*autoclavable

## Digital implant analogue \*\*

Item code
● WDA29SI
● WDA33SI



## Digital analogue for Performing Abutment

Item code
DAPRA



\*\*Libraries can be downloaded from [www.globalwin.eu](http://www.globalwin.eu)

# Prosthetic abutments

Abutments made from:

**titanium, sterilise before use**


VP2SI screw for  
prosthetic abut-  
ment



## Winbase abutments \*\*

for single-unit screw-retained CAD CAM prosthesis, VP2SI included


Gingival height	Item code
GH 1.0 mm	● WB2901SI
GH 2.0 mm	● WB2902SI
GH 3.0 mm	● WB2903SI
GH 1.0 mm	● WB3301SI
GH 2.0 mm	● WB3302SI
GH 3.0 mm	● WB3303SI



## Premilled Abutments\*\*

SLIM premilled abutments VP2SI prosthetic screw included  
Slim Internal Hex Connection Attachment M


Diameter	Item code
Ø 10.00 mm	WPM10MSI
Ø 12.00 mm	WPM12MSI




\*\*Libraries can be downloaded from [www.globalwin.eu](http://www.globalwin.eu)



# GLOBALWIN Guided surgery





**BIOSAFIN**<sup>®</sup>  
**CREA** DIGITAL CENTER

### **SPECIALISED DIGITAL CENTRE**

Guided surgery is one of the services offered by the innovative CREA Digital Centre that BIOSAFIN has set up within its own production site in Trezzano Rosa (Milan).

The CREA Digital Centre is a state-of-the-art specialist centre dedicated to those wishing to work with a partner that puts research and development, expertise, support, quality and precision at the forefront to offer tailored solutions and enhance the clinician's skills and the patient's restorations.

With significant investments in production technology, consolidated expertise and by always keeping abreast of new developments, the CREA Digital Centre is able to provide the following services:

- GUIDED SURGERY
- CUSTOMISATION OF COMPONENTS AND PROSTHETIC STRUCTURES
- MANUFACTURE OF BIALIGNER INVISIBLE ALIGNERS

### **FULLY CONTROLLED MANUFACTURE**

BIOSAFIN manufacturing is certified:

- ISO 9001: customer-satisfaction oriented
- ISO 13485: specifically relevant to the quality of medical devices

The innovative aspect and the quality of the devices have gained considerable international recognition and clinical testimonies, which supports the sound scientific background on which the company bases its long history of reliability.


# Guided surgery today

## Simple and precise procedure planning.

The inclusion of Guided Surgery protocols is ever more a part of dental surgical practice as these help optimise procedure times and results for restorations on implants, with great patient satisfaction as well.

Once the diagnostic examination has been performed, this approach makes it possible to proceed immediately to the planning of the restoration, guided by software that provides all of the necessary tools.

The intervention can be performed using the flapless surgical technique – without the need to raise a mucoperiosteal flap – and thus proceed directly to placement of the temporary restoration.



## APPROVED PROTOCOLS

By planning the procedure in 3D with CREA - 3D Guided Surgery software, it is possible to consider various restoration solutions, simulate prosthetic solutions and thoroughly examine the individual functional and aesthetic parameters. Restoration times are thus optimised and surgical risks are minimised, with results that are in line with expectations.

## GUIDED SURGERY IMPLANT PROSTHESIS RESTORATION WORKFLOW



## GUIDED SURGERY KIT INSTRUMENTS

Item Code	Description
STD0G	Mucotome for guided surgery
STD16G	6mm-long Pilot drill for guided surgery
DSD12	12mm-long Short manual screwdriver for torque wrench
DSD16	16mm-long Medium manual screwdriver for torque wrench
HSD20	20mm-long Contra-angle screwdriver
HMG	Contra-angle adaptor for Guided Surgery with MTRG series instruments
ADEXG	Extension for Guided Surgery for use with the ADMTRG adaptor
STD28G	Drill for guided surgery Ø: 2.2/2.6 mm - L: 8 mm
STD210G	Drill for guided surgery Ø: 2.2/2.6 mm - L: 10 mm
STD2115G	Drill for guided surgery Ø: 2.2/2.6 mm - L: 11.5 mm
STD213G	Drill for guided surgery Ø: 2.2/2.6 mm - L: 13 mm
STD215G	Drill for guided surgery Ø: 2.2/2.6 mm - L: 15 mm
STD38G	Drill for guided surgery Ø: 2.6/3.0 mm - L: 8 mm
STD310G	Drill for guided surgery Ø: 2.6/3.0 mm - L: 10 mm
STD3115G	Drill for guided surgery Ø: 2.6/3.0 mm - L: 11.5 mm
STD313G	Drill for guided surgery Ø: 2.6/3.0 mm - L: 13 mm
STD315G	Drill for guided surgery Ø: 2.6/3.0 mm - L: 15 mm
STD48G	Drill for guided surgery Ø: 3.0/3.4 mm - L: 8 mm
STD410G	Drill for guided surgery Ø: 3.0/3.4 mm - L: 10 mm
STD4115G	Drill for guided surgery Ø: 3.0/3.4 mm - L: 11.5 mm
STD413G	Drill for guided surgery Ø: 3.0/3.4 mm - L: 13 mm
STD415G	Drill for guided surgery Ø: 3.0/3.4 mm - L: 15 mm
STD58G	Drill for guided surgery Ø: 3.4/3.8 mm - L: 8 mm
STD510G	Drill for guided surgery Ø: 3.4/3.8 mm - L: 10 mm
STD5115G	Drill for guided surgery Ø: 3.4/3.8 mm - L: 11.5 mm
STD513G	Drill for guided surgery Ø: 3.4/3.8 mm - L: 13 mm
STD515G	Drill for guided surgery Ø: 3.4/3.8 mm - L: 15 mm
STD68G	Drill for guided surgery Ø: 3.8/4.2 mm - L: 8 mm
STD610G	Drill for guided surgery Ø: 3.8/4.2 mm - L: 10 mm
STD6115G	Drill for Guided surgery Ø: 3.8/4.2 mm - L: 11.5 mm
STD613G	Drill for guided surgery Ø: 3.8/4.2 mm - L: 13 mm
STD615G	Drill for guided surgery Ø: 3.8/4.2 mm - L: 15 mm
BP35G	Guided Cortical Drill Ø 3.40 mm
BP37G	Guided Cortical Drill Ø 3.65 mm
BP43G	Guided Cortical Drill Ø 4.20 mm
BP50G	Guided Cortical Drill Ø 4.90 mm
STDPG	1.5 mm-long Pin drill for guided surgery
PGG	GLOBALWIN Surgical stent fastening pin
DTW75	Torque wrench 10-75 Ncm
ADMTRG	Adaptor for guided surgery for DTW75 torque wrench
MMTRG	Mounter for guided surgery
GWGST	GLOBALWIN® Surgical tray for guided surgery

## CBCT ANATOMICAL LANDMARKS

SGG30	Landmark spheres for Guided Surgery (3.0 mm) Pack of 20
US003-01	Universal stent pack of 1 (CBCT anatomical landmarks)
US003-03	Universal stent pack of 3 (CBCT anatomical landmarks)
US003-05	Universal stent pack of 5 (CBCT anatomical landmarks)
US003-10	Universal stent pack of 10 (CBCT anatomical landmarks)

## SURGICAL STENTS

GG501I	1 implant Male Surgical Stent Lower - Bushings included
GG501S	1 implant Male Surgical Stent Upper - Bushings included
GG502I	2 implants Male Surgical Stent Lower - Bushings included
GG502S	2 implants Male Surgical Stent Upper - Bushings included
GG503I	3 implants Male Surgical Stent Lower - Bushings included
GG503S	3 implants Male Surgical Stent Upper - Bushings included
GG504I	4 implants Male Surgical Stent Lower - Bushings included
GG504S	4 implants Male Surgical Stent Upper - Bushings included
GG505I	5 implants Male Surgical Stent Lower - Bushings included
GG505S	5 implants Male Surgical Stent Upper - Bushings included
GG506I	6 implants Male Surgical Stent Lower - Bushings included
GG506S	6 implants Male Surgical Stent Upper - Bushings included
GG507I	7 implants Male Surgical Stent Lower - Bushings included
GG507S	7 implants Male Surgical Stent Upper - Bushings included
GG508I	8 implants Male Surgical Stent Lower - Bushings included
GG508S	8 implants Male Surgical Stent Upper - Bushings included
GG517I	Additional Surgical Stent - Lower - Bushings included
GG517S	Additional Surgical Stent - Upper - Bushings included

## SOFTWARE AND SERVICES

M03013DML-CRE-V	CREA-3D Guided surgery Viewer
M03013DML-CRE	CREA-3D Guided surgery Software
C3D0SC	One Shot CREA-3D - 30 min consultation included - Virtual Patient creation included
C3DOS	One Shot CREA-3D Virtual Patient creation included
C3DVP	CREA-3D Virtual Patient creation

## ACCESSORIES

Item Code	Description
BGG50	GLOBALWIN Guided Sleeve 5.0 mm
HMTRSG	SLIM hex connection mounter for guided surgery



**Guided Surgery Kit**



# GLOBALWIN

IMPLANT SYSTEM

GLOBALWIN medical devices are compliant with EC Directive 93/42 as amended



# Rehabilitative technique for multi-unit screw-retained prostheses


## Teeth Just on 4

In order to offer satisfactory and economically viable restorative solutions to a significant portion of the population, we have developed specific surgical and prosthetic protocols for resolving cases with extensive edentulism the upper or lower jaws:

### Teeth Just on 4 and Teeth Just on 6.

Clinical studies and research have shown that the inclined insertion of two distal implants is an effective and simple technique to compensate for any bone failure, thus expanding the prosthetic base in a stable and functional manner.

In the case of the mandible, this technique also protects the final portion of the arch, avoiding interference with the mandibular nerve. The Teeth Just On 4 or Teeth Just on 6 technique allows complete, permanent and stable arch rehabilitation, in many cases avoiding the necessity of bone regeneration procedures and the discomfort and cost to the patient involved.




**TEETH JUST ON 4 - GLOBALWIN Implants**  
For the complete, permanent and stable prosthetic rehabilitation of the lower arch on only 4 implants.

## CAB Clip Abutment Bar

The **CAB** CLIP Abutment BAR device promotes the solidifying of dental implants by means of a Clip - Abutment - Bar combination applied to immediate and non-immediate loading, screw-retained prostheses, in line with the Teeth Just on 4/6 restorative technique protocol.

The CAB is used to build a passive titanium structure, to be implemented in the shortest time, in the context of immediate load application.

It represents a reinforcement of the temporary prosthesis, minimising the risk of fracturing that could lead to failure of the Implants.



USA

The patent issued for the CAB® by the US authority confirms its originality and innovative technological content.



European Union

International and European Patent  
PCT/EP2011/072448  
EP Patent no. 11425032.7

## Bar preparation procedure

After placing the implants and their respective abutments from the PRAxxxx line, an impression is taken with the appropriate PRAA abutments.

Subsequently, the laboratory model is developed with analogues for PRAxxxx fitted. The AT... abutments are selected according to the shoulder height, allowing the bar to be parallel to the occlusal plane. Once the ATxx abutments are in place, each individual bar is cut using the appropriate CB instrument.

After cutting all the bars, the appropriate CF1 or CM1 Clip is inserted, and the entire structure is assembled by mounting it on the ATxx abutments and permanently fixing it to the latter with appropriate cements.

## CAB Kit




## CAB Components


Code BTO Bar without Clip	
Code BT1 Retentive bar with Male clip, cylindrical hole and 1mm thickness	
Code BT2 Retentive bar with Male clip, cylindrical hole and 2mm thickness	
Code BT3 Retentive bar with male clip with elliptical hole and 2 mm thickness	
Code CM1 Male Clip with cylindrical hole and 1 mm thickness	
Code CM2 Male Clip with cylindrical hole and 2 mm thickness	
Code CM3 Male Clip with elliptical hole and 2 mm thickness	
Code CF1 Female Clip	
Code AT. CAB abutment shoulder where xx indicates variable shoulder height from 1.7 to 4.2 mm	
Code CB Cutter Bar	

## Fully adjustable bar

Thanks to the cylindrical or elliptical geometry of the special clips, the Clip-Abutment connection allows extreme versatility of use even in the event of severe disparallelism. This feature makes it ideal for multi-unit screw-retained prostheses with immediate or deferred loading.



1 Measurement of the CAB bar with the help of the Cutter Bar device



2 Cutting the CAB bar with separating disc



3 The mounted CAB bar is easily positioned, even in the event of severe disparallelism

# CAB Kit composition

To carry out a first CAB case, it is necessary to have a complete CAB Kit.

After the first case, there will be a minimum of 20 to a maximum of 50 Abutments available in the KIT.

For subsequent cases, it is sufficient to reorder only the missing items.

All kits include ATxx Abutments for each shoulder height from 1.7 mm to 4.2 mm

Just On 4 KIT		
Code	Description	Quantity
CM1	Male Clip – Cylindrical Hole 1 mm Thickness	1
CF1	Female Clip – Cylindrical Hole	2
BT1	Bar + Cylindrical CM 1 mm thickness	1
BT2	Bar + Cylindrical CM 2 mm thickness	2
BT3	Bar + Elliptical CM 2 mm thickness	2
AT17	CAB Abutment - Shoulder h 1.7mm	4
AT22	CAB Abutment - Shoulder h 2.2mm	4
AT27	CAB Abutment - Shoulder h 2.7mm	4
AT32	CAB Abutment - Shoulder h 3.2mm	4
AT37	CAB Abutment - Shoulder h 3.7mm	4
AT42	CAB Abutment - Shoulder h 4.2mm	4
PRAA	PRAxxxx line analogue	4
CB	CAB Cutter Bar	2
TCAB	Abutment Bar Clip Tray	1
DSD16	16mm-long Medium manual screwdriver for torque wrench	1

Just On 6 KIT		
Code	Description	Quantity
CM1	Male Clip – Cylindrical Hole 1 mm Thickness	1
CF1	Female Clip – Cylindrical Hole	4
BT1	Bar + Cylindrical CM 1 mm thickness	3
BT2	Bar + Cylindrical CM 2 mm thickness	2
BT3	Bar + Elliptical CM 2 mm thickness	2
AT17	CAB Abutment - Shoulder h 1.7mm	6
AT22	CAB Abutment - Shoulder h 2.2mm	6
AT27	CAB Abutment - Shoulder h 2.7mm	6
AT32	CAB Abutment - Shoulder h 3.2mm	6
AT37	CAB Abutment - Shoulder h 3.7mm	6
AT42	CAB Abutment - Shoulder h 4.2mm	6
PRAA	PRAxxxx line analogue	6
CB	CAB Cutter Bar	2
TCAB	Abutment Bar Clip Tray	1
DSD16	16mm-long Medium manual screwdriver for torque wrench	1

Just On 4/6 KIT		
Code	Description	Quantity
CM1	Male Clip – Cylindrical Hole 1 mm Thickness	2
CF1	Female Clip – Cylindrical Hole	6
BT1	Bar + Cylindrical CM 1 mm thickness	4
BT2	Bar + Cylindrical CM 2 mm thickness	4
BT3	Bar + Elliptical CM 2 mm thickness	4
AT17	CAB Abutment - Shoulder h 1.7mm	10
AT22	CAB Abutment - Shoulder h 2.2mm	10
AT27	CAB Abutment - Shoulder h 2.7mm	10
AT32	CAB Abutment - Shoulder h 3.2mm	10
AT37	CAB Abutment - Shoulder h 3.7mm	10
AT42	CAB Abutment - Shoulder h 4.2mm	10
PRAA	PRAxxxx line analogue	10
CB	CAB Cutter Bar	2
TCAB	Abutment Bar Clip Tray	1
DSD16	16mm-long Medium manual screwdriver for torque wrench	1

Items available to purchase separately	
Code	Description
CM1	Male Clip – Cylindrical Hole 1 mm Thickness
CM2	Male Clip – Cylindrical Hole 2 mm Thickness
CM3	Male Clip – Elliptical Hole 2 mm Thickness
CF1	Female Clip – Cylindrical Hole
BT0	Bar without Clip L. 30 mm - divisible
BT1	Bar + Cylindrical CM 1 mm thickness
BT2	Bar + Cylindrical CM 2 mm thickness
BT3	Bar + Elliptical CM 2 mm thickness
AT17	CAB Abutment - Shoulder h 1.7mm
AT22	CAB Abutment - Shoulder h 2.2mm
AT27	CAB Abutment - Shoulder h 2.7mm
AT32	CAB Abutment - Shoulder h 3.2mm
AT37	CAB Abutment - Shoulder h 3.7mm
AT42	CAB Abutment - Shoulder h 4.2mm
PRAA	PRAxxxx line analogue
CB	CAB Cutter Bar
TCAB	Abutment Bar Clip Tray
DSD16	16mm-long Medium manual screwdriver for torque wrench

# Materials

## Metals

Type	Composition
Titanium	Ti grade 4 ASTM F67
Titanium	ASTM F132 Grade 5 Ti
Stainless Steel	AISI 303 - AISI 630 - ASTM F899
Platinum Alloy	Au 60%, Pd 15%, Pt 24,9%, Ir 0,1% Melting range 1350°-1460° C
CrCo	Cobalt chrome molybdenum alloy compliant to standard ISO5832-12: 2019/ASTM F1537

## Synthetic Resins

Type	Composition
Polymethylmethacrylate	PMMA Melting range 150–160°C

# Sterilisation protocol



- Detergent used: SEKUSEPT or similar, dilution 1 measure of detergent per litre of water.
- Decontaminate in solution for 15 min.
- Ultrasonic cleaning for 15 min. at 60°C.
- Dry (very important).

**NOTE:** After checking for the presence of any residues, the equipment must be bagged and sterilised at 134°C 2 ATM - 2 min. - (Europe); 135°C/270°F - 3 min. - (USA).

*Equipment processed in an autoclave must be rinsed and dried thoroughly since the autoclaving process increases the oxidising action of detergents.*

# Patient Information:



## Implant Passport

An important document reporting procedures undergone and materials used, thanks to their batch number on an adhesive sticker: an additional guarantee of quality and safety, should it be necessary in the future to identify any of them.

# Select bibliography:

## **Implant prosthetic rehabilitation in HIV-positive patients: a comparison of two different implant surface roughness**

M. C. Francia, S. Galli, F. Bova, S. Ferrari Parabita, P. Capparè, E. F. Gherlone  
Journal of Osseointegration 2017; 9(1): 77

## **Combined microcomputer tomography, biomechanical and histomorphometric analysis of the peri-implant bone: a pilot study in minipig model**

M. Gramanzini, S. Gargiulo, F. Zarone, R. Megna, A. Apicella, R. Aversa, M. Salvatore, M. Mancini, R. Sorrentino, A. Brunetti  
Dental Materials 32 (2016) 794-806

## **NEWSLETTER BIOSAF IN: JUST ON 4: 4 ANNI DOPO: GENNAIO 2011 – APRIL 2015**

Authors: Operational Unit of Dentistry - Università Vita - Salute San Raffaele of Milan, Dir. Prof. Enrico Gherlone

## **Edentulous patients rehabilitated according to the "all-on-four" procedure with prefabricated bar system**

P. Capparè, M. Pasi, G. L. Di Domenico, L. Cisternino, E. Polizzi, R. Crespi - Department of Dentistry, IRCCS San Raffaele Hospital, Milan, Italy- Dental School, Vita-Salute University, Milan, Italy Minerva Stomatologica 2015; 64 (Suppl.1 al No2): 81

## **Crestal bone remodeling around platform switched, immediately loaded implants placed in sites of previous failures**

A. Quaranta, A. Cicconetti, L. Battaglia, M. Piemontese, G. Pompa, I. Vozza  
European Journal of Inflammation Vol.10 N°2 Accepted July 30 2012

## **Isolation of osteogenic progenitors from amniotic fluid using a single-step cell culture protocol**

S. Tetè, I. Antonucci, G. D'Apolito, U. Di Tore, V. Zizzari, A. De Carlo, M. Tumedei, L. Stuppia, E. Gherlone  
Doctor Os, supplemento maggio 2011 - XXII (5)

## **"Isolation of osteogenic progenitors from human amniotic fluid using a single step culture protocol"**

I. Antonucci, I. Iezzi, E. Morizio, F. Mastrangelo, A. Pantalone, M. Mattioli Belmonte, A. Gigante, V. Salini, G. Calabrese, S. Tetè, G. Palka, L. Stuppia  
SILENCE a Journal of RNA regulation, 2009

## **"In vitro behaviour onto different titanium surface of osteoblast-like cells obtained from human dental pulp"**

S. Tetè, F. Mastrangelo, V. Zizzari, G. D'Apolito, N. Fiorentino, U. Desiato, M.T. Sberna, R. Quaresima, L. Stuppia, R. Vinci, E. F. Gherlone  
Atti del 7th Annual Meeting of ISSCR Int. Society of Stem Cell Research, Barcellona July 2009

## **"Novel Protocol of osteogenic differentiation from amniotic fluid cells"**

S. Tetè, F. Mastrangelo, M. Tranasi, V. Zizzari, I. Antonucci, G. D'Apolito, T. Marchese, R. Vinci, L. Stuppia, E. F. Gherlone  
Atti del 7th Annual Meeting of ISSCR Int. Society of Stem Cell Research, Barcellona July 2009

## **"Interfaccia osso - impianto nei differenti tipi di carico degli impianti dentali"**

S. Tetè, G. D'Apolito, F. Mastrangelo, R. Vinci, E. F. Gherlone  
Atti III Expo di Autunno Università della Lombardia, 27-28 Novembre 2009

## **"Valutazione della capacità osteogenica di hafscs ottenute da liquido amniotico"**

S. Tetè, U. Di Tore, V. Zizzari, L. Stuppia, F. Zarone, E. F. Gherlone  
Atti III Expo di Autunno Università della Lombardia, 27-28 Novembre 2009

## **"In vitro evaluation of osteoblast-like cells from different sources"**

F. Mastrangelo, M. Tranasi, V. Zizzari, D. Farronato, T. Traini, R. F. Grassi, L. Stuppia, S. Tetè  
Atti del 86th General Session of International Association of Dental Research (I.A.D.R.) Toronto, 2-5 Giugno 2008

## **"A macro-and nanostructure evaluation of a Novel Dental Implant"**

S. Tetè, F. Mastrangelo, T. Traini, R. Vinci, G. Sammartino, G. Marenzi, E. F. Gherlone  
Implant Dentistry vol.17n.3 2008

## **"Valutazione del comportamento di osteoblasti derivanti da cellule staminali di liquido amniotico su differenti superfici implantari"**

A. Desiderio, A. D'Incecco, L. Stuppia, I. Antonucci, S. Tetè  
Atti I Expo di Autunno Università della Lombardia, 1 Dicembre 2007

## **"Evaluation of surface nano-topographic effect on afscells growth"**

A. Desiderio, F. Mastrangelo, T. Traini, L. Stuppia, G. Sammartino, S. Tetè  
Atti dell'Annual Meeting del I.A.D.R. International Association of Dental Research, Continental Division Thessaloniki 26-29 Settembre 2007

## **Impianto troncoconico con componente protesica sottodimensionata per mantenere l'osso creatale**

L. Prosper, S. Redaelli, A. D'Addona, E. F. Gherlone  
Poster 13° Congresso Nazionale del Collegio dei Docenti di Odontoiatria - Roma 5-8 Aprile 2006

## **Four-year follow-up of larger diameter implants placed in fresh extraction sockets using a resorbable membrane or a resorbable alloplastic material**

JOMI The L. Prosper, E. F. Gherlone, S. Redaelli, M. Quaranta  
JOMI The International Journal of Oral & Maxillofacial Implants 2003; vol.18 - n.6

## **Analysis of 435 screw-vent dental implants placed in 161 patients: software enhancement of clinical evaluation**

Artin ML.  
Implant Dent. 2002; 11 (1): 58-66

## **"Osteointegrazione a 5 anni di impianti sabbiati o sabbiati e mordenzati ritenenti protesi parziali fisse"**

L. Prosper, F. Di Carlo, I. Vozza, M. Quaranta  
QI Quintessence International 1 - 2002

## **"Clinical trial on osseointegration using sandblasted or sandblasted and acidetched implants"**

S. Redaelli, L. Prosper, F. Di Carlo, A. Daddona, M. Quaranta - 2001

## **"Healing period of titanium implants with sandblasted and acid-etched surface"**

L. Prosper, F. Di Carlo, S. Redaelli, R. Scaringi, M. Quaranta  
3° World Congress of Osteointegration 2001

## **"Tempo di guarigione ossea per impianti di titanio con superficie sabbiata o superficie sabbiata e mordenzata"**

L. Prosper, F. Di Carlo, S. Redaelli, G. Radaelli, R. Scaringi, M. Quaranta  
La rivista Internazionale di Odontoiatria Protetica, vol.13 n. 1 - 2000

## **"Analisi in vivo su tre superfici implantari: valutazione istomorfometrica"**

S. Redaelli, L. Prosper, F. Di Carlo, A. Scarano, I. Vozza - anno 2000

BIOSAFIN

Partner of



Making Your Life Better.

**Manufactured by:**

**BIOSAF IN** srl  
info@biosafin.com  
www.biosafin.com

**Offices in Italy:**

**MILAN:**

Via Cagliari 32/44  
20060 – Trezzano Rosa – (MI)  
Zona Industriale  
Tel. +39 02 90968692  
Fax +39 02 90968541

**ANCONA:**

Via Tiraboschi, 36/G  
60131 - Ancona (AN)  
Tel. +39 071 2071897  
Fax +39 071 203261

CODE 19 1 / ITA / VI - 2024



**GLOBALWIN**  
I M P L A N T   S Y S T E M

GLOBALWIN medical devices are compliant with EC Directive 93/42 as amended

[www.globalwin.eu](http://www.globalwin.eu)



**BIOSAF IN** is certified to:

**ISO 9001**, which certifies the entire work process, from start to finish, further demonstrating our compliance with quality standards considered optimal for the protection of product Users – Professionals – and end users – Patients.

**ISO 13485** specifically relevant to the quality of Medical Devices.